
# Baylis

MED TECH

## Passport™

*Steerable Sheath*

en- Instructions for Use • Steerable Sheath

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fr- Mode d'emploi • Gaine orientable

This page is intentionally left blank.

# **Baylis** MED TECH

# **Passport**<sup>TM</sup> *Steerable Sheath*


en- English..... Page 2

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fr- Français..... Page 18

## Explanation of symbols on package labeling

Refer to the outer package label to see which symbols apply to this product.

	Consult Instructions for use	<b>Rx only</b>	By prescription only
	Instructions for Use	<b>REF</b>	Catalogue Number
	Do not reuse	<b>LOT</b>	Lot number
	Do not resterilize!		Length
	Do not use if package is damaged!		Inside diameter
<b>STERILE EO</b>	Sterilized using ethylene oxide		Outer diameter
	Single sterile barrier system	<b>MAX GW</b>	Maximum Recommended Guidewire
	Keep away from sunlight		Curve
	Keep Dry		Package Contains
	Use-by date	<b>QTY</b>	Quantity
	Date of manufacture		Distributor
	Manufacturer		Uni-Directional Curve
	Caution		

# English • Table of Contents

1. Device Description.....	3
2. Intended Use/Indications for Use.....	3
3. Cautions.....	3
4. Contraindications.....	4
5. Adverse Effects.....	4
6. Possible Complications.....	4
7. Precautions.....	4
8. Necessary Hospital Equipment.....	5
9. Directions for Use.....	5
10. Storage.....	8
11. Handling.....	8
12. Disposal.....	8
13. U.S. Limited Warranty and Disclaimer.....	8

---

## 1. Device Description

The Passport™ Steerable Sheath is a steerable guiding sheath that features an adjustable tip through use of a rotating collar that allows the physician to deflect the distal tip section of the sheath. The sheath has a hemostatic valve, sideport with stopcock, radiopaque distal tip, contains side flush portholes (available in some models), and the tip is deflected by rotating the collar incorporated into the handle. The dilator has depth markings (on some models) and snap or screw locks onto the sheath (depending on the model). The dilator has a tapered distal tip and an inner lumen.

### Contents of Package (Please refer to product label for contents):

1. Steerable Guiding Sheath
2. Dilator
3. Guidewire (optional)
4. Valve Bypass Tool (optional)

## 2. Intended Use/Indications for Use

The Passport™ Steerable Sheath is intended for the introduction of diagnostic and therapeutic devices into the human vasculature, including but not limited to intracardiac, renal or other peripheral placements. Do not use this device for neural placements.

## 3. Cautions

- Device is supplied sterile. Do not use if package has been previously opened or damaged.
- Prior to use, read all package inserts, warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.
- Dimensional specifications of the device are located on product labels.
- Procedure must be performed by trained medical personnel well versed in anatomical landmarks, safe technique, and potential complications.
- The product is designed for single use only. Risks of reusing the product are: Structural integrity might get compromised and/or biological tissue residue and other microorganisms that may infect a patient.
- Do not re-sterilize or reuse. Do not alter the device in any way.

- Do not use this device for neural placements.
- If a serious incident has occurred and it is suspected to be attributed to the Passport™ Steerable Sheath, please contact Oscor Inc. at 727-937-2511 or email us at rga@Oscor.com.

#### 4. Contraindications

- Known active systemic or local infection
- Known inability to obtain vascular access
- Patients with atrial thrombus or myxoma, or interatrial baffle or patch
- Use of a steerable sheath is contraindicated in patients with obstructed or inadequate vasculature.

#### 5. Adverse Effects

Potential adverse effects related to the use of the steerable sheath include but are not limited to:

- |                                       |   |
|---------------------------------------|---|
| • Air Embolism                        | • Intimal tear                              |
| • Allergic reaction to contrast media | • Irregular heart beat                      |
| • Aortic puncture                     | • Local nerve damage                        |
| • Arrhythmias                         | • Mediastinal widening                      |
| • Arteriovenous fistula formation     | • Myocardial infarction                     |
| • Atrial septal defect                | • Pacemaker/defibrillator lead displacement |
| • Bleeding plexus injury              | • Perforation                               |
| • Catheter entrapment                 | • Pericardial/pleural effusion              |
| • Cardiac tamponade                   | • Pneumothorax                              |
| • Coronary artery spasm and/or damage | • Pseudoaneurysm formation                  |
| • Vessel spasm                        | • Pulmonary edema                           |
| • Dislodgement                        | • Stroke                                    |
| • Dissection                          | • Subclavian artery puncture                |
| • Endocarditis                        | • Thromboembolic events                     |
| • Heart Block                         | • Thrombophlebitis                          |
| • Hematoma formation                  | • Valve damage                              |
| • Hemorrhage                          | • Vascular occlusion                        |
| • Hemothorax                          | • Vasovagal reaction                        |
| • Infection                           | • Vessel damage/Vessel trauma               |
|                                       | • Death                                     |

#### 6. Possible Complications

Arrhythmias may occur during the use of any intracardiac device. Careful monitoring and availability of emergency equipment are mandatory.

#### 7. Precautions

**Transvenous device compatibility:** Use the steerable sheath only with compatible transvenous devices. Use the appropriate size sheath for the size of the transvenous device being utilized. Consequences of using the steerable sheath with incompatible devices may include the inability to deliver the transvenous device or damage to the transvenous device during delivery.

#### Before you flush the sheath:

- Do not connect a power injection syringe to the sideport and inject contrast solution.
- Do not aspirate steerable sheath with a guidewire in place through the hemostatic valve.
- Aspiration with a guidewire through the valve may cause air embolism which can result in significant morbidity or death.
- Use the sideport, located on the handle, to inject contrast solution or flush the steerable sheath. The steerable sheath must be thoroughly flushed with either saline or heparinized saline and free of air prior to use to avoid air embolism to the patient.
- Aspiration and flushing of the sheath should be performed frequently to help minimize the potential for air embolism.
- For injecting or aspirating through the sheath, use the sideport only with stopcock. Prior to infusion, remove all air using the sideport.
- Indwelling sheath should be internally supported by a catheter, electrode, or dilator.

- Never advance, torque, or withdraw sheath when resistance is met. Determine cause by fluoroscopy and then take remedial action.
- The following conditions require special care if a transseptal approach is used:
  - » enlarged aortic root
  - » marked right atrial enlargement
  - » small left atrium
  - » marked distortion of the thoracic configuration (for example, kyphosis or scoliosis)

## 8. Necessary Hospital Equipment

- Keep external defibrillation equipment nearby for immediate use during insertion, placement, acute lead system testing, or whenever arrhythmias are possible or intentionally induced.
- If the patient has left bundle branch block, back up pacing should be readily available during insertion of the steerable sheath assembly. Use of steerable sheath assembly may cause heart block.

## 9. Directions for Use

**Note:** The steerable sheath is designed for temporary use (contemplated implant duration of eight hours or less) within the peripheral, renal, coronary and intracardiac systems. Proper surgical procedures and sterile techniques are the responsibility of the medical professional. The following procedures are provided for information only. Each physician must apply the information in these instructions according to professional medical training and experience.

The procedure for using the steerable sheath to gain peripheral, renal, coronary and intracardiac system access include the following steps:

- 9.1 Preparing the steerable sheath for insertion
- 9.2 General use of the steerable sheath
- 9.3 Transseptal use of the steerable sheath
- 9.4 Deflecting and straightening the steerable sheath
- 9.5 Removing the steerable sheath

### 9.1 Preparing Steerable Sheath for Insertion:

**Warning:** The steerable sheath and dilator must be thoroughly flushed with either saline or heparinized saline and free of air prior to use to avoid air embolism to the patient.

1. Remove the steerable sheath and sheath dilator from the sterile package within the sterile field.
2. Verify deflecting and straightening of the distal section of the steerable sheath using the handle of the sheath. Refer to section 9.4 “Deflecting and straightening the steerable sheath”.
3. Thoroughly flush steerable sheath and dilator with either saline or heparinized saline. The steerable sheath and dilator must be free of air prior to use to avoid air embolism.

### 9.2 General Use of the Steerable Sheath:

**Warning:** If the patient has left bundle branch block, back up pacing should be readily available during insertion of the steerable sheath assembly. Use of the steerable sheath assembly may cause heart block.

1. Select appropriate size of the sheath for the device to be delivered.
2. Follow standard access procedures to the desired branch within the vasculature.
3. Wet the dilator shaft with sterile saline solution prior to insertion through the hemostatic valve. Note: Any diagnostic device, therapeutic device or component inserted through the hemostatic valve of the sheath should be wet and placed through the center of the valve to prevent tearing of the seal and leakage.
4. When appropriate, insert the guidewire according to its own Instructions for Use.
5. Assemble the dilator and sheath together and snap or screw on the dilator hub to the sheath hub prior to insertion.

6. Insert the sheath assembly over the guidewire to access the desired branch.
7. Do not force the steerable sheath assembly if significant resistance is encountered during the insertion or passage. If resistance is encountered, determine the cause and correct before continuing the procedure.
8. Aspirate all air from the sheath by connecting the syringe to the sideport with stopcock.
9. When in desired position, unsnap or unscrew the dilator from the sheath and slowly pull it back to prevent air ingress.
10. Refer to section 9.4 “Deflecting and straightening the steerable sheath” for instructions.
11. Remove the dilator slowly when access in desired vasculature is reached.
12. Remove the guidewire slowly.
13. Deliver the diagnostic and/or therapeutic device through the sheath hemostatic valve following its own Instructions for Use.
14. The sheath sideport with stopcock should be used throughout the procedure to aspirate or deliver liquids.

**NOTE:** Any diagnostic device, therapeutic device or component inserted through the hemostatic valve of the sheath should be wet and placed through the center of the valve to prevent tearing of the seal and leakage. In some cases, a valve bypass tool (VBT) might be used to ease the insertion of the diagnostic and/or therapeutic device through the valve. Improper use of a valve bypass tool (VBT) can cause air embolism and back bleeding.

### **9.3 For Transseptal Use of the Steerable Sheath:**

**Warning:** If the patient has left bundle branch block, back up pacing should be readily available during insertion of the steerable sheath assembly. Use of the steerable sheath assembly may cause heart block.

1. Follow the Instructions for Use that were supplied with the transseptal needle (not supplied) and the guidewire. If a transseptal needle is used, it is recommended to use together with a stylet in order to prevent skiving of the inner lumen of dilator.
2. Select appropriate size of the sheath for the device to be delivered and appropriate sheath length to match the transseptal needle. Refer to transseptal needle Instructions for Use for needed sheath length.
3. Wet the dilator shaft with sterile saline solution prior to insertion through the hemostatic valve.
4. Place dilator inside the sheath and snap or screw on both hubs together. Thread the dilator/sheath assembly over the guidewire, using a slight twisting motion. Note: Any device/component inserted through the hemostatic valve of the sheath should be wet and placed through the center of the valve to prevent tearing of the seal and leakage.
5. Aspirate all air from the sheath by connecting the syringe to the sideport.
6. Advance the sheath/dilator assembly into the superior vena cava, just above the right atrium.
7. Separate the dilator and sheath hub approximately 1 cm. This will aid the introduction of the curved transseptal needle.
8. Slowly remove the guidewire from the dilator. Remove all air from the dilator by slowly aspirating blood and then flush the dilator. Then completely flush the transseptal needle.
9. Introduce the properly prepared transseptal needle into the dilator hub. Carefully advance the curved section of the needle into the dilator, being certain not to restrict the movement of the needle.
10. Re-attach the dilator and sheath hubs.
11. While maintaining the sheath position, slowly advance the curved section of the needle until it is about to protrude from the dilator tip. Good right atrial pressure should be observed before proceeding.

12. Position the needle and sheath set in the right atrium. Verify position under fluoroscopy.
13. Position dilator and needle point against the atrial septum in the region of the fossa ovalis by gradually rotating the needle posterior and toward the left scapula during withdrawal. Use continual pressure monitoring and repeated anterior-posterior and lateral visualization of the tip under fluoroscopy during all positioning procedures.
14. Advance the needle and complete the transseptal puncture. Measure the left atrium pressure immediately after needle penetration. Do not advance the dilator if the pressure is not acceptable.
15. Advance the dilator with the needle in place through the septum. Acceptable left atrial pressure should be continually observed. A sequential increase in resistance to movement followed by a sharp decrease in resistance will indicate location of the dilator in the left atrium.
16. Withdraw the needle point with the dilator tip in place. The dilator with the needle point within should be freely located in the left atrium. Verify with fluoroscopy.
17. Advance the sheath slowly over the dilator-needle combination until it is in the left atrium. A slow rotating motion of the sheath as firm pressure is applied will aid in this procedure. The sheath will be in position when a sharp reduction of resistance is felt.
18. Advance the sheath approximately 2 cm into the left atrium while maintaining the dilator/needle position.
19. After disconnecting the pressure monitoring line from the needle, slowly remove the needle from the dilator.
20. Slowly remove the dilator from the sheath. Caution: Rapid removal may damage the valve membrane resulting in blood flow through the valve.
21. Attach the sheath sideport to the monitoring line. Gently aspirate blood through the sideport for sampling and to make sure the sheath is clear of air. Blood should aspirate freely through the sideport.
22. Under fluoroscopy, frequently observe the location of the radiopaque tip marker to maintain the location of the sheath in the left atrium.
23. Refer to section 9.4 “Deflecting and Straightening the Steerable Sheath” for instructions.
24. Introduce the catheter through the hemostatic valve into the left atrium. Improved catheter manipulation may be obtained by withdrawing the sheath into the right atrium. The sheath should be returned to the left atrium over the catheter prior to removing the catheter.
25. For removal of the sheath please refer to section 9.5 “Removing of the Steerable Sheath”.
26. After removal of the sheath, use standard technique to achieve hemostasis.

#### **9.4 Deflecting and Straightening the Steerable Sheath:**

The Passport™ Steerable Sheath features a deflecting distal tip for site specific placement of the sheath and subsequent device(s) within the peripheral, renal, coronary and intracardiac systems. The deflection mechanism is controlled by operator manipulation of the rotating collar handle.

1. Twist the deflection collar slowly and carefully to deflect the distal tip. **Caution:** The sheath will deflect at the speed which the deflection collar is turned. Avoid rapid deflection that may cause vessel damage.
2. When the desired deflection point or site access is achieved, stop turning the deflection collar.

**Caution:** To ensure retention of the desired deflection position(s), avoid contact with the deflection collar that may cause the sheath to change deflection shape. Do not deflect the sheath with dilator or with the device inserted.

 **Caution:** When in the deflected position, the steerable guiding sheath should not be rotated while encountering resistance since that could severely damage the vessel, heart chamber, or anatomy of the patient.

### **9.5 Removing the Steerable Sheath:**

 **Warning:** Before removing the steerable sheath, straighten the distal section to avoid damage during removal. Refer to section 9.4 “Deflecting and Straightening the Steerable Sheath” for instructions.

1. Remove any device from the steerable sheath prior to retracting the sheath.
2. Straighten the distal tip by rotating the collar to the straight position.
3. Slowly retract the device from the body, if resistance is met, stop retracting the device and ensure that the tip is in the straight position.
4. Discard the steerable sheath after it has been removed from the body in accordance with institution’s biohazard requirements.

## **10. Storage**

Store product in a cool, dark and dry place. Do not expose product to excessive light, including sun light or ultraviolet light, or excessive humidity.

## **11. Handling**

- Avoid subjecting the device to unusual stresses. Handle the steerable sheath with care at all times.
- Before removing the steerable sheath from packaging, straighten the distal section as much as possible to avoid damage during removal. Refer to section 9.4 “Deflecting and straightening the steerable sheath” for instructions.
- Do not kink, stretch, or severely bend the steerable sheath.
- Do not use surgical instruments to grasp the steerable sheath.
- Do not use excessive force when inserting the steerable sheath into a vessel.
- Do not force the steerable sheath if significant resistance is encountered during insertion or passage.
- Do not wipe the sheath with chlorofluorocarbon based solvents.
- Ensure the steerable sheath is thoroughly flushed and free of air prior to use.
- Avoid contact with liquids other than isopropyl alcohol, blood, saline or contrast solution.
- An unsupported sheath (without an inserted device or dilator) may be susceptible to kinking during advancement or torsional manipulation.

## **12. Disposal**

Discard of all used devices in accordance with institution’s biohazard requirements.

## **13. U.S. Limited Warranty and Disclaimer**

Oscor Inc. hereby warrants that if any Oscor Inc. product fails to perform within normal tolerances due to a defect in materials or workmanship, Oscor Inc. will provide, at no charge, a replacement of Oscor Inc. product. This limited warranty applies only if each of the following conditions are met:

- a. The product was designed and manufactured by Oscor Inc.
- b. The product was stored according to label information.
- c. The product was used according to its intended use.
- d. The product was used by trained and certified clinical staff.
- e. The product has not been mishandled, reprocessed, or altered in any way.
- f. The product was used before the “Use-by date” marked on the packaging of the product (when applicable).


No representation or warranty is made that an Oscor Inc. product will not fail. Oscor Inc. disclaims responsibility for any medical complication, including death, resulting from the use of its products. Except as expressly provided by this limited warranty, Oscor Inc. is not responsible for any direct, incidental, or consequential damages based on any defect, failure,

or malfunction of its products whether the claim is based on warranty, contract, tort, or otherwise. Some states do not allow the exclusion or limitation of consequential damages and so the above limitation or exclusion may not apply to Buyer.

Except as expressly provided by the limited warranty, Oscan Inc. makes no warranty, express or implied, including but not limited to, any implied warranty of merchantability, or fitness for a particular purpose.

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Passport™ and the Baylis Medtech logo are trademarks and/or registered trademarks of Baylis Medical Technologies Inc. in the United States of America and/or other countries.
- a. Le produit a été conçu et fabriqué par Oscor Inc.
## **Baylis Medical Technologies Inc.**
Palm Harbor, FL 34683 U.S.A.