Instructions For Use - Nykanen Radiofrequency Wire

Content Effective Date: Jun 15, 2023

Content Revision: 20

Originator: Jonathan Jebaraj

PRODUCT

Nykanen Radiofrequency Wire

Instructions for Use Revision Information

| Document Number | Revision | Change Request Number |
|-----------------|----------|-----------------------|
| DMR-00691 | 19 | 05820 |
| DMR-00691 | 20 | 07474 |

Begin Instructions for Use on following page

Printed Date: 18-Jul-2023

Page 1 of 26

| | | | | | |
|----------------------|-------------------------------------|--------------------|-----|-------------------|---|
| TEMPLATE INFORMATION | Template Document Number: QMS-00415 | Template CR Number | N/A | Template Revision | 1 |
|----------------------|-------------------------------------|--------------------|-----|-------------------|---|

CONFIDENTIAL © 2021-2023: The contents of this document are confidential and are intended solely for use by Baylis Medical Technologies personnel. This document must not be copied, used or disclosed to a third party without the written consent of Baylis Medical Technologies. Document is only a controlled document if stamped with "Doc Control Released" and date. All other printed documents are only valid on date printed.

Instructions for Use

Nykanen Radiofrequency Wire

| | |
|-----------------|-----|
| [English | 1] |
| [Français..... | 3] |
| [Deutsch..... | 5] |
| [Italiano..... | 7] |
| [Español..... | 9] |
| [Português..... | 12] |
| [Čeština..... | 14] |
| [Dutch..... | 16] |
| [Dansk..... | 18] |
| [Suomi..... | 20] |
| [Norsk..... | 22] |
| [Svenska..... | 24] |

BAYLIS MEDICAL TECHNOLOGIES INC.
2645 Matheson Blvd. East
Mississauga, ON, Canada L4W 5S4
+1 (905) 948-5800
www.baylismedtech.com

© Copyright [Baylis Medical Technologies Inc., 2023]

The Baylis MedTech logo
is a trademark and/or registered trademark of Baylis Medical Technologies Inc.
in the United States of America and/or other countries.

English

Carefully read all instructions prior to use. Observe all warnings and precautions noted in these instructions. Failure to do so may result in patient complications. Baylis Medical Technologies relies on the physician to determine, assess and communicate to each individual patient all foreseeable risks of the procedure.

I. DEVICE DESCRIPTION

The Nykanen Radiofrequency Wire delivers radiofrequency (RF) power in a monopolar mode between its distal electrode and a commercially available external Disposable Indifferent (Dispersive) Patch (DIP) Electrode, which is in compliance with IEC 60601-2-2:2017 :2017. The Nykanen Wire is connected to the BMC Radiofrequency Puncture Generator via the Catheter Connector Cable. The cable can accommodate a maximum allowable length of 3 meters. Detailed information concerning the BMC Radiofrequency Puncture Generator is contained in a separate manual that accompanies the Generator (entitled "BMC Radiofrequency Puncture Generator Instructions for Use"). BMC Generators compatible with the Nykanen Radiofrequency Wire include the RFP-100 or RFP-100A.

The critical dimensions of the Nykanen Radiofrequency Wire can be found on the device label. The Teflon insulation on the wire body facilitates smooth advancement through guiding catheters, while the outer diameter is small enough to provide adequate spacing for injection of contrast solution through the guiding catheter. The distal portion of the Nykanen Radiofrequency Wire is flexible and the active tip is shaped to be atraumatic to the cardiac tissue unless RF energy is applied.

II. INDICATIONS FOR USE

U.S.A.: The Nykanen Radiofrequency Wire is intended to cut and/or coagulate soft tissues.

Canada: The Nykanen Radiofrequency Wire is intended to create an atrial septal defect in the heart.

III. CONTRAINDICATIONS

U.S.A.: The Nykanen Radiofrequency Wire is not recommended for use with any conditions that do not require cutting or coagulation of soft tissues.

Canada: The Nykanen Radiofrequency Wire is not recommended for use with any conditions that do not require the creation of an atrial septal defect.

IV. WARNINGS

- The Nykanen Radiofrequency Wire is supplied STERILE using an ethylene oxide process. Do not use if the package is damaged.
- Laboratory staff and patients can undergo significant x-ray exposure during radiofrequency puncture procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure.
- The Nykanen Radiofrequency Wire is intended for single patient use only. Do not attempt to sterilize and reuse the wire. Reuse can cause the patient injury and/or the communication of infectious disease(s) from one patient to another. Failure to do so may result in patient complications.
- Do not attempt to use the Nykanen Radiofrequency Wire with electrocautery devices.
- The Nykanen Radiofrequency Wire must be used with the Catheter Connector Cable and the BMC Radiofrequency Puncture Generator. Attempts to use it with other RF Generators and connector wires can result in electrocution of the patient and/or operator.
- The active tip of the Nykanen Radiofrequency Wire is fragile. Be careful not to damage the tip while handling the wire. If the tip becomes damaged discard the wire immediately.
- If the active tip of the wire becomes bent at any time during its use, dispose of the wire immediately. Do not attempt to straighten the active tip.

- For RFP-100: Do not attempt to puncture with an initial power setting of greater than 10 watts. The initial attempt should be made with a 10 watts setting. In subsequent attempts, the power setting can be increased, if necessary.
- During RF output, implanted devices such as pacemakers may be affected. Qualified advice should be obtained as necessary, to minimize the risk from injury due to implanted device malfunction.
- The Nykanen Radiofrequency Wire is not intended for use with neonatal patients (less than one month of age). Do not attempt to treat neonatal patients with the wire.
- Individual patient anatomy and physician technique may require procedural variations.
- Inadvertent activation or movement of the activated electrode outside of the field of vision may result in injury to the patient.

V. PRECAUTIONS

- Do not attempt to use the Nykanen Radiofrequency Wire or ancillary equipment before thoroughly reading the accompanying Instructions for Use.
- Radiofrequency puncture procedures should be performed only by physicians thoroughly trained in the techniques of radiofrequency powered puncture in a fully equipped catheterization laboratory.
- The sterile packaging should be visually inspected prior to use to detect any compromise. Ensure that the packaging has not been damaged. Do not use the equipment if the packaging has been compromised.
- Visually inspect the Nykanen Radiofrequency Wire prior to use to ensure there is no cracking or damage to the insulating material. Do not use the wire if there is any damage.
- Do not use the Nykanen Radiofrequency Wire after the "Use By" date indicated on the label.
- The Nykanen Radiofrequency Wire is intended for use with only those devices listed in section VIII, "Equipment Required."
- Read and follow the manufacturer's instructions for use of the Disposable Indifferent (Dispersive) Patch (DIP) electrode. Always use DIP electrodes that meet or exceed IEC 60601-2-2:2017 requirements.
- Placement of the dispersive electrode on the thigh could be associated with higher impedance.
- In order to prevent the risk of ignition make sure that flammable material is not present in the room during RF power application.
- Take precautions to limit the effects that the electromagnetic interference (EMI) produced by the Generator may have on the performance of other equipment. Check the compatibility and safety of combinations of other physiological monitoring and electrical apparatus to be used on the patient in addition to the Generator.
- Adequate filtering must be used to allow continuous monitoring of the surface electrocardiogram (ECG) during radiofrequency power applications.
- Careful wire manipulation must be performed to avoid cardiac damage, or tamponade. Wire advancement should be done under fluoroscopic guidance. If resistance is encountered, DO NOT use excessive force to advance or withdraw the wire.
- Do not attempt to puncture until firm position of the active tip has been achieved against the atrial septum.
- Do not activate RF while tip is within a third party introducer, dilator or catheter.
- It is recommended not to exceed five (5) radiofrequency power applications per Nykanen Radiofrequency Wire.
- Do not bend the Nykanen Radiofrequency Wire. Excessive bending or kinking of the wire shaft may damage the integrity of the wire and may cause patient injury. Care must be taken when handling the wire.
- The RF Generator is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the wire and DIP electrode, particularly when operating the device.
- During power delivery, the patient should not be allowed to come in contact with ground metal surfaces.
- Apparent low power output or failure of the equipment to function properly at normal settings may indicate faulty application of the DIP electrode, failure to an electrical lead, or poor tissue contact at the active tip. Check for obvious equipment defects or misapplication. Attempt to better position the tip of the Nykanen Radiofrequency Wire against the atrial septum. Only increase the power if low power output persists.
- Regularly inspect and test re-usable cables and accessories.
- Baylis Medical Technologies relies on the physician to determine, assess and communicate to each individual patient all foreseeable risks of the Baylis Medical Radiofrequency Puncture System.
- Do not attempt to insert and use the proximal end of the Nykanen Radiofrequency Wire as the active tip.
- The Nykanen Radiofrequency Wire and Catheter Connector Cable are rated for a maximum voltage of 400 rms.

VI. ADVERSE EVENTS

Adverse events which may occur while using the Baylis Medical Radiofrequency Puncture System include:

| | |
|-------------------------------|--------------------------------------|
| Atrial Fibrillation (AF) | Ventricular Tachycardia (VT) |
| Atrial Flutter | Vessel Spasm |
| Atrial Tachycardia | Perforation and/or Cardiac Tamponade |
| Sustained Arrhythmia | Sepsis/Infection/Inflammation |
| Embolic Event | Pericardial/Pleural Effusion |
| Vessel Dissection/Puncture | Pain and Tenderness |
| Hemorrhage | Conduction/Heart Block |
| Hematoma | Pericarditis |
| Ventricular Fibrillation (VF) | Death |
| Cerebral Infarction | Myocardial Infarction |

VII. INSPECTION PRIOR TO USE


Prior to use of the Baylis Medical Radiofrequency Puncture System, the individual components including the BMC Radiofrequency Puncture Generator, Nykanen Radiofrequency Wire, and the Catheter Connector Cable should be carefully examined for damage or defects, as should all equipment used in the procedure. Inspect and ensure no exposed metal along the body of the device and between the proximal end and the push-button connector of the catheter connector cable. Inspect all connections before use. Do not use defective equipment. Do not reuse the Nykanen Radiofrequency Wire.

VIII. EQUIPMENT REQUIRED

Intracardiac puncture procedures should be performed in a specialized clinical setting equipped with a fluoroscopy unit, radiographic table, physiologic recorder, emergency equipment and instrumentation for gaining vascular access. Ancillary materials required to perform cardiac puncture include:

- BMC Radiofrequency Puncture Generator
- Catheter Connector Cable (Model RFP-101 for use with RFP-100 Generator or RFX-BAY-OTW-10-SU for use with RFP-100A Generator)

- Disposable Indifferent (dispersive) Patch (DIP) electrode that meets or exceeds IEC 60601-2-2:2017 requirements for electrosurgical electrodes, such as the Valley Labs Polyhesive Electrode #E7506.
 - It is recommended to use Nykanen Radiofrequency Wire with Baylis Medical ProTrack™ Microcatheter (Model CIC 38-145).
- IX. SETTING UP THE SYSTEM**
- Please refer to the BMC Radiofrequency Puncture Generator Instructions for Use. The Instructions for Use describe the steps to follow to connect the system, set parameters, and deliver radiofrequency power.
 - All instructions for equipment required should be carefully read, understood, and followed. Failure to do so may result in complications.
- X. ATTACHING THE INDIFFERENT ELECTRODE**
- Remove the Disposable Indifferent (dispersive) Patch (DIP) electrode from the packaging and peel off the backing to expose the conductive gel surface. Check to be sure the pad is moist and sticky to the touch before placing it on the patient. A dry electrode will have limited grounding capability.
 - Place the DIP electrode on a well-vascularized convex skin surface which is in close proximity to the heart. Do not place this electrode on the thigh, since this location is associated with higher impedance (see section V, "Precautions"). Avoid scar tissue, body prominence, adipose tissue, and areas where fluid may pool. Shave, clean, and dry the application site as needed.
 - Check to be sure that excellent contact has been achieved over the entire area of the DIP electrode. Burns can result when RF power is delivered to a DIP electrode with poor contact.
 - Plug the DIP electrode connector into the appropriate port on the BMC RF Puncture Generator as per the Generator Instructions for Use. Make sure that the DIP electrode connector is firmly pressed into the socket.
- XI. SUGGESTED DIRECTIONS FOR USE**
- The Nykanen Radiofrequency Wire is usually inserted into the femoral vein and is then positioned into the appropriate chamber of the heart under fluoroscopic guidance. Once appropriate positioning has been achieved, radiofrequency power is delivered via the BMC RF Puncture Generator. This results in the puncture of the targeted cardiac tissue.

- XII. CLEANING AND STERILIZATION INSTRUCTIONS**
The Nykanen Radiofrequency Wire is single use and supplied sterile. Do not clean or re-sterilize the device. The devices can be considered sterile only if the package is not opened or damaged prior to use.
- XIII. DISPOSAL OF WASTE**
Treat used Nykanen Radiofrequency Wire as biohazard and dispose of in compliance with standard hospital procedures.
- XIV. CUSTOMER SERVICE AND PRODUCT RETURN INFORMATION**
If you have any problems with or questions about Baylis Medical Equipment contact our technical support personnel.

BAYLIS MEDICAL TECHNOLOGIES INC.
2645 Matheson Blvd. East
Mississauga, ON, Canada L4W 5S4
+1 (905) 948-5800
www.baylismedtech.com

- NOTES:**
1. In order to return products you must have a return authorization number before shipping the products back to Baylis Medical Technologies. Product Return Instructions will be provided to you at this time.
 2. Ensure that any product being returned to Baylis Medical has been cleaned, decontaminated and/or sterilized as indicated in the Product Return Instruction before returning it for warranted service. Baylis Medical will not accept any piece of used equipment that has not been properly cleaned or decontaminated as per the Product Return Instructions.

XV. LABELING AND SYMBOLS

| | |
|----------------|---|
| | Manufacturer |
| | Sterile using ethylene oxide |
| | Use By |
| | Caution |
| | Follow Instructions for Use |
| | Model number |
| Rx ONLY | Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician. |
| | Single Use – Do not reuse |
| | Do Not Resterilize |
| | Lot Number |
| | Do Not Use if Packaging is Damaged |
| | Keep Away From Sunlight |
| | Non-Pyrogenic |
| | Use with BMC RFP Generators Only |
| | Maximum of five (5) RF Power Applications |

- XVI. LIMITED WARRANTIES**
Baylis Medical Technologies Inc. (BMT) warrants its Disposable and Accessory products against defects in materials and workmanship. BMT warrants that sterile products will remain sterile for a period of time as shown on the label as long as the original package remains intact. Under this

Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMT will replace or repair, in its absolute and sole discretion, any such product, less any charges to BMT for transportation and labor costs incidental to inspection, removal or restocking of product. The length of the warranty is: (i) for the Disposable products, the shelf life of the product, and (ii) for the Accessory products, 90 days from shipment date. This limited warranty applies only to new original factory delivered products that have been used for their normal and intended uses. BMT's Limited Warranty shall not apply to BMT products which have been resterilized, repaired, altered, or modified in any way and shall not apply to BMT products which have been improperly stored or improperly cleaned, installed, operated or maintained contrary to BMT's instructions.

DISCLAIMER AND EXCLUSION OF OTHER WARRANTIES

Baylis Medical Technologies Inc. (BMT) warrants its Disposable and Accessory products against defects in materials and workmanship. BMT warrants that sterile products will remain sterile for a period of time as shown on the label as long as the original package remains intact. Under this Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMT will replace or repair, in its absolute and sole discretion, any such product, less any charges to BMT for transportation and labor costs incidental to inspection, removal or restocking of product. The length of the warranty is: (i) for the Disposable products, the shelf life of the product, and (ii) for the Accessory products, 90 days from shipment date. This limited warranty applies only to new original factory delivered products that have been used for their normal and intended uses. BMT's Limited Warranty shall not apply to BMT products which have been resterilized, repaired, altered, or modified in any way and shall not apply to BMT products which have been improperly stored or improperly cleaned, installed, operated or maintained contrary to BMT's instructions.

LIMITATION OF LIABILITY FOR DAMAGES

THE LIMITED WARRANTY ABOVE IS THE SOLE WARRANTY PROVIDED BY SELLER. SELLER DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE OR PURPOSE.

THE REMEDY SET FORTH HEREIN SHALL BE THE EXCLUSIVE REMEDY FOR ANY WARRANTY CLAIM, AND ADDITIONAL DAMAGES, INCLUDING CONSEQUENTIAL DAMAGES OR DAMAGES FOR BUSINESS INTERRUPTION OR LOSS OF PROFIT, REVENUE, MATERIALS, ANTICIPATED SAVINGS, DATA, CONTRACT, GOODWILL OR THE LIKE (WHETHER DIRECT OR INDIRECT IN NATURE) OR FOR ANY OTHER FORM OF INCIDENTAL, OR INDIRECT DAMAGES OF ANY KIND, SHALL NOT BE AVAILABLE. SELLER'S MAXIMUM CUMULATIVE LIABILITY RELATIVE TO ALL OTHER CLAIMS AND LIABILITIES, INCLUDING OBLIGATIONS UNDER ANY INDEMNITY, WHETHER OR NOT INSURED, WILL NOT EXCEED THE COST OF THE PRODUCT(S) GIVING RISE TO THE CLAIM OR LIABILITY. SELLER DISCLAIMS ALL LIABILITY RELATIVE TO GRATUITOUS INFORMATION OR ASSISTANCE PROVIDED BY, BUT NOT REQUIRED OF SELLER HEREUNDER. ANY ACTION AGAINST SELLER MUST BE BROUGHT WITHIN EIGHTEEN (18) MONTHS AFTER THE CAUSE OF ACTION ACCRUES. THESE DISCLAIMERS AND LIMITATIONS OF LIABILITY WILL APPLY REGARDLESS OF ANY OTHER CONTRARY PROVISION HEREOF AND REGARDLESS OF THE FORM OF ACTION, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND FURTHER WILL EXTEND TO THE BENEFIT OF SELLER'S VENDORS, APPOINTED DISTRIBUTORS AND OTHER AUTHORIZED RESELLERS AS THIRD-PARTY BENEFICIARIES. EACH PROVISION HEREOF WHICH PROVIDES FOR A LIMITATION OF LIABILITY, DISCLAIMER OF WARRANTY OR CONDITION OR EXCLUSION OF DAMAGES IS SEVERABLE AND INDEPENDENT OF ANY OTHER PROVISION AND IS TO BE ENFORCED AS SUCH.

IN ANY CLAIM OR LAWSUIT FOR DAMAGES ARISING FROM ALLEGED BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, PRODUCT LIABILITY OR ANY OTHER LEGAL OR EQUITABLE THEORY, THE BUYER SPECIFICALLY AGREES THAT BMT SHALL NOT BE LIABLE FOR DAMAGES OR FOR LOSS OF PROFITS, WHETHER FROM BUYER OR BUYER'S CUSTOMERS. BMT'S LIABILITY SHALL BE LIMITED TO THE PURCHASE COST TO BUYER OF THE SPECIFIED GOODS SOLD BY BMT TO BUYER WHICH GIVE RISE TO THE CLAIM FOR LIABILITY.

No agent, employee or representative of Baylis Medical has the authority to bind the Company to any other warranty, affirmation or representation concerning the product.

This warranty is valid only to the original purchaser of Baylis Medical products directly from a Baylis Medical authorized agent. The original purchaser cannot transfer the warranty.

Use of any BMT product shall be deemed acceptance of the terms and conditions herein.

The warranty periods for Baylis Medical products are as follows:

| | |
|---------------------|--------------------------------|
| Disposable Products | The shelf life of the product |
| Accessory Products | 90 days from the shipment date |

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