The PowerWire® Pulse Radiofrequency Guidewire Kit includes the PowerWire® Pulse Radiofrequency Guidewire and the RFP-100A Connector Cable. Carefully read all instructions prior to use.

Instructions for Use

PowerWire® Pulse Radiofrequency Guidewire Kit

English	1
Espanol	3

Baylis Medical Technologies Inc.
2645 Matheson Blvd. East
Mississauga, Ontario, Canada, L4W 5S4
+1 (905) 948-5800
www.baylismedtech.com

© Baylis Medical Company Inc. 2007-2021. Amendments © Baylis Medical Technologies Inc. 2023.

PowerWire and the Baylis MedTech logo
are trademarks and/or registered trademarks of Baylis Medical Technologies Inc.

English

Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these instructions. Failure to do so may result in patient complications.

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

I. DEVICE DESCRIPTION

The PowerWire® Pulse Radiofrequency Guidewire Kit includes the PowerWire® Pulse Radiofrequency Guidewire and the RFP-100A Connector Cable.

The PowerWire® Pulse Radiofrequency Guidewire must be used with an approved Baylis Radiofrequency Puncture Generator (Baylis RF Generator) and Baylis Connector Cable.

The PowerWire® Pulse Radiofrequency Guidewire delivers radiofrequency (RF) power in a monopolar mode between its distal electrode and a commercially available external Disposable Indifferent (Dispersive) Patch (DIP) Electrode, which is in compliance with the IEC 60601-2-2 standards. The PowerWire® Pulse Radiofrequency Guidewire is connected to the Baylis RF Generator via the Baylis Connector Cable.

Dimensions for different models of the PowerWire® Pulse Radiofrequency Guidewire are described on the product label. The device consists of an insulated core nitinol wire with a rounded tip. The insulation on the PowerWire® Pulse Radiofrequency Guidewire body provides electrical insulation and lubricity. The distal portion of the PowerWire® Pulse Radiofrequency Guidewire is flexible and the active tip is rounded to be atraumatic to vessel walls.

The RFP-100A Connector Cable connects the RFP-100A Baylis Radiofrequency Puncture Generator (RFP-100A Generator) to the PowerWire® Pulse Radiofrequency Guidewire. This Cable enables radiofrequency (RF) power to be delivered from the Generator to a PowerWire® Pulse Radiofrequency Guidewire.

Detailed information concerning the RFP-100A Generator is contained in a separate manual that accompanies the Generator (RFP-100A Generator Instructions for Use).

The dimensions for the PowerWire® Pulse Radiofrequency Guidewire and the RFP-100A Connector Cable can be found on the device label and in Section VI "Product Specifications." The RFP-100A Connector Cable has a multi-pin connector on one end that mates with the RFP-100A Generator and a connector at the other end, which mates with a PowerWire® Pulse Radiofrequency Guidewire.

II. INDICATIONS FOR USE

The PowerWire® Pulse Radiofrequency Guidewire Kit is indicated for use to cut and/or coagulate soft tissue.

III. CONTRAINDICATIONS

The PowerWire® Pulse Radiofrequency Guidewire Kit is not recommended for use with any conditions that do not require cutting and/or coagulation of soft tissue.

The RFP-100A Connector Cable included in the PowerWire® Pulse Radiofrequency Guidewire Kit is not recommended for use with any other RF generator or any other device.

IV. WARNINGS

- Only physicians with a thorough understanding of angiography and percutaneous interventional procedures should use the PowerWire® Pulse Radiofrequency Guidewire Kit. It is recommended that physicians avail themselves of pre-clinical training, a review of pertinent literature and other appropriate education before attempting new interventional procedures.
- The PowerWire® Pulse Radiofrequency Guidewire is supplied STERILE using an ethylene oxide process. Do not use if the package is damaged.
- The PowerWire® Pulse Radiofrequency Guidewire and RFP-100A Connector Cable are intended for single patient use only. Do not attempt to sterilize and reuse either device. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another.
- The PowerWire RF Guidewire must be used with the Baylis Connector Cable. Attempts to use it with other connector cables can result in electrocution of the patient and/or operator.
- The RFP-100A Connector Cable must only be used with the RFP-100A Generator and the PowerWire® Pulse RF Guidewire. Attempts to use it with other RF Generators and devices can result in electrocution of the patient and/or operator.
- To perform electrosurgery safely and effectively, the surgeon and operating room staff must be properly trained in the surgical techniques that are to be employed.
- Do not use electrosurgery in the presence of flammable anesthetics or other flammable gases, near flammable fluids or objects, or in the presence of oxidizing agents, as fire could result.
- Place active accessories in a clean, dry, non-conductive, and highly visible area away from the patient when not in use. Inadvertent contact with the patient may result in burns. Contact with drapes may cause fire.

- The active tip of the PowerWire® Pulse Radiofrequency Guidewire is fragile. Be careful not to damage the tip while handling the PowerWire® Pulse RF Guidewire. If the tip becomes damaged, discard the PowerWire® Pulse RF Guidewire immediately.
- If the active tip of the PowerWire® Pulse Radiofrequency Guidewire becomes bent at any time during its use, dispose of the PowerWire® Pulse Radiofrequency Guidewire immediately. Do not modify the PowerWire® Pulse Radiofrequency Guidewire tip. Modifications to the tip may result in tip breakage or other damage.
- Laboratory staff and patients can undergo significant x-ray exposure during radiofrequency cutting procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure.

V. PRECAUTIONS

- Do not attempt to use the PowerWire® Pulse Radiofrequency Guidewire Kit before thoroughly reading the accompanying Instructions for Use.
- Radiofrequency interventional procedures should be performed only by physicians thoroughly trained in the technique in an operating room fully equipped for interventional procedures and/or a fully equipped catheterization laboratory.
- The sterile packaging and PowerWire® Pulse Radiofrequency Guidewire Kit should be visually inspected prior to use to detect any compromise. Ensure that the packaging has not been damaged. Check the PowerWire® Pulse Radiofrequency Guidewire and the RFP-100A Connector Cable to ensure there is no cracking or damage to the insulating material.
- Do not use the PowerWire® Pulse Radiofrequency Guidewire Kit after the "Use By" date indicated on the label.
- Do not bend the PowerWire® Pulse Radiofrequency Guidewire and the RFP-100A Connector Cable. Excessive bending or kinking of the PowerWire® Pulse Radiofrequency Guidewire shaft may damage the integrity of the device and may cause patient injury. Care must be taken when handling the PowerWire® Pulse Radiofrequency Guidewire and the RFP-100A Connector Cable.
- Take precautions to limit the effects that the electromagnetic interference (EMI) produced by the Generator may have on the performance of other equipment. Check the compatibility and safety of combinations of other physiological monitoring and electrical apparatus to be used on the patient in addition to the Generator.
- Adequate filtering must be used to allow continuous monitoring of the surface electrocardiogram (ECG) during radiofrequency power applications.
- Ensure that the Disposable Indifferent (Dispersive) Patch (DIP) electrode is appropriately selected, properly applied to the patient, and connected to the Baylis RF Generator. Read and follow the manufacturer's instructions for use of the DIP electrode. Always use DIP electrodes that meet or exceed current IEC 60601-2-2 requirements.
- Placement of the dispersive electrode on the thigh could be associated with higher impedance.
- The Generator is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the PowerWire® Pulse Radiofrequency Guidewire and/or DIP electrode, particularly when operating the device.
- During power delivery, the patient should not be allowed to come in contact with ground metal surfaces.
- In order to prevent the risk of ignition, make sure that flammable material is not present in the room during RF power application.
- Confirm that the PowerWire® Pulse Radiofrequency Guidewire is properly connected to the correct power output receptacle on the RFP-100A Generator.
- Examine all devices to be connected to the RFP-100A Generator. After connection, ensure that they are functioning as intended.
- Confirm proper RFP-100A Generator settings before proceeding with surgery.
- Careful manipulation of the PowerWire® Pulse Radiofrequency Guidewire must be performed to avoid tissue trauma. PowerWire® Pulse Radiofrequency Guidewire advancement should be done under fluoroscopic guidance. If resistance is encountered, DO NOT use excessive force to advance or withdraw the PowerWire® Pulse Radiofrequency Guidewire.
- Do not attempt to deliver radiofrequency energy until the tip of the PowerWire® Pulse Radiofrequency Guidewire is confirmed to be in good contact with the target biological material.
- It is recommended not to exceed twenty (20) radiofrequency power applications per PowerWire® Pulse Radiofrequency Guidewire.
- During power delivery, the patient should not be allowed to come in contact with ground metal surfaces.
- Care must be taken to verify the compatibility of the PowerWire® Pulse Radiofrequency Guidewire when used with catheters.
- Do not activate RF whilst tip is within catheter.
- The active tip may remain hot enough to cause burns after the electrosurgical current is deactivated.

- Inadvertent activation or movement of the activated electrode outside of the field of vision may result in injury to the patient.
- Localized burns to the patient or physician may result from electrical currents carried through conductive objects. Electrical current may be generated in conductive objects by direct contact with the active electrode being in close proximity to the conductive object.
- Never disconnect the RFP-100A Connector Cable from the RFP-100A Generator while the Generator is delivering RF power.
- Never disconnect the RFP-100A Connector Cable from the RFP-100A Generator by pulling on the cable. Failure to disconnect the cable properly may result in damage to the cable.
- Do not twist the RFP-100A Connector Cable while inserting or removing it from the Isolated Patient Connector on the RFP-100A Generator. Twisting the cable may result in damage to the pin connectors.
- Apparent low power output or failure of the equipment to function properly at normal settings may indicate faulty application of the DIP electrode, failure to an electrical lead, or poor contact with the target biological material at the active tip. Check for obvious equipment defects or misapplication. Attempt to better position the active tip of the PowerWire® Pulse Radiofrequency Guidewire against the target biological material.
- Baylis Medical Technologies relies on the physician to determine, assess and communicate to each individual patient all foreseeable risks of the PowerWire® Pulse Radiofrequency Guidewire.

VI. PRODUCT SPECIFICATIONS

Product	PowerWire® Pulse Radiofrequency Guidewire	Product	RFP-100A Connector Cable
Useable Length	250cm	Useable Length	10 feet (3m)
Outside Diameter	0.035" / 0.89 mm	Generator Connector	multi-pin
Marker Bands	5 Placed 1 cm apart	Device Connector	Push Button

The PowerWire RF Guidewire is a defibrillation-proof type CF applied part.

VII. ADVERSE EVENTS

Adverse events that may occur while performing the radiofrequency cutting procedure include:

Thrombosis	Pain and Tenderness
Vessel perforation	Distal embolization
Vessel dissection	Sepsis/Infection
Hemorrhage	Hematoma
Effusion	Arrhythmias
Tamponade	Death

VIII. EQUIPMENT REQUIRED

Radiofrequency cutting procedures should be performed in a specialized clinical setting equipped with a fluoroscopy unit, radiographic table, physiologic recorder, emergency equipment and instrumentation for gaining vascular access. Ancillary materials required to perform this procedure include:

- Baylis Radiofrequency Generator
- Disposable Indifferent (dispersive) Patch (DIP) electrode

IX. INSPECTION PRIOR TO USE

Prior to performing the radiofrequency cutting procedure, the individual components including the Baylis Radiofrequency Generator, PowerWire® Pulse Radiofrequency Guidewire, and the RFP-100A Connector Cable should be carefully examined for damage or defects, as should all equipment used in the procedure.

Inspect the packaging of the PowerWire® Pulse Radiofrequency Guidewire Kit to ensure the package has not been damaged, sterility has not been compromised, and that a product label is present. Ensure that the RFP-100A Connector Cable and the PowerWire® Pulse Radiofrequency Guidewire have no visible damage, such as discoloration, cracks, label fading, splicing, or kinks. Do not use damaged or defective equipment. Do not reuse the PowerWire® Pulse Radiofrequency Guidewire or the RFP-100A Connector Cable.

X. DIRECTIONS FOR USE

All instructions for equipment required should be carefully read, understood, and followed. Failure to do so may result in complications.

- The PowerWire® Pulse Radiofrequency Guidewire Kit is supplied sterile. Use aseptic technique when opening the package and handling the product in the sterile field.
- Connect the generator connector end of the cable to the isolated patient connector port on the RFP-100A Generator as per the Generator Instructions for Use. The RFP-100A Connector Cable uses a circular connector, keyed for proper alignment. Gently line up the connector pins with the socket and

push in until the connector fits firmly into the socket. Any attempt to connect the cable otherwise will damage the pins on the connector.

- Do not use excessive force in connecting the cable to the generator. Use of excessive force may result in damage to the connector pins.
- Firmly grasp the catheter connector end of the cable in one hand. Using your thumb depress the red button on the top of the connector. Slowly insert the proximal end of the RF puncture wire into the opening of the catheter connector. Once the exposed portion of the proximal end of the device is no longer visible release the red button on the connector. Gently tug on the device to ensure that you have a secure connection.
- Once the desired site has been accessed and the tip of the PowerWire® Pulse Radiofrequency Guidewire is in contact against the tissue, the PowerWire® Pulse Radiofrequency Guidewire can be used to therapeutically deliver RF power. RF power is delivered to the tissue via the distal tip electrode, which results in thermal necrosis (cutting and coagulation) of the tissue. To deliver RF power from the generator to the PowerWire® Pulse Radiofrequency Guidewire, follow the instructions provided with the generator.
- Apply firm pressure to the PowerWire® Pulse Radiofrequency Guidewire during the application of radiofrequency energy to successfully advance it through the target site.
- Radiofrequency power delivery can be terminated by pressing the RF ON/OFF button on the Generator if the timer has not expired.
- Successful advancement of the PowerWire® Pulse Radiofrequency Guidewire can be confirmed by monitoring the PowerWire® Pulse Radiofrequency Guidewire under fluoroscopy.
- If advancement through the target site is not successful after twenty (20) radiofrequency power applications, it is advised that the user utilize an alternate method for the procedure.
- To disconnect the PowerWire® Pulse Radiofrequency Guidewire from the RFP-100A Connector Cable, depress the red button on the catheter connector and gently remove the proximal end of the puncture device from the connector cable.
- To disconnect the cable from the generator, grasp the connector firmly and gently pull it straight out of the socket.

XI. CLEANING AND STERILIZATION INSTRUCTIONS

The PowerWire® Pulse Radiofrequency Guidewire Kit is intended for single use only. Do not clean or re-sterilize any device in the PowerWire® Pulse Radiofrequency Guidewire Kit. Do not use solvents to clean the distal tip of the PowerWire® Pulse Radiofrequency Guidewire.

The PowerWire® Pulse RF Guidewire Kit can be considered sterile only if the package is not opened or damaged prior to use.

XII. CUSTOMER SERVICE AND PRODUCT RETURN INFORMATION

If you have any problems with or questions about Baylis Medical Technologies equipment, contact technical support personnel.

NOTES:

- 1 In order to return products, you must have a return authorization number before shipping the products back to Baylis Medical Technologies.
- 2 Baylis Medical Technologies will not accept any piece of used equipment without a sterilization certificate. Ensure that any product being returned to Baylis Medical Technologies has been cleaned, decontaminated and sterilized as per user instructions before returning it for warranted service.

XIII. TROUBLESHOOTING


The following table is provided to assist the user in diagnosing potential problems.

PROBLEM	COMMENTS	TROUBLESHOOTING
Connector Cable does not fit into the Isolated Patient Connector on the front panel of the generator	The connectors are designed to connect in a specific way for safety reasons. If the connector "keys" are out of line, the connectors won't fit together	Check that the connector keys are lined up in the proper orientation. Ensure that the connectors are clean and unobstructed.
Generator Error Messages	In order to successfully cut biological material using RF energy, all devices must be properly connected and in good working order.	Ensure that all connections are made: - PowerWire® Pulse Radiofrequency Guidewire to Connector Cable - Connector Cable to Generator - Generator to power outlet - Generator to DIP electrode - DIP electrode to patient Visually inspect the PowerWire® Pulse Radiofrequency Guidewire

		or Connector Cable for damage. Immediately discard any damaged equipment. If problem persists discontinue use. For error messages encountered while attempting radiofrequency cutting, refer to the Instructions for Use document that accompanies the Generator. If errors persist, attach a new Connector Cable. If this solves the problem, discard the damaged connector cable.
PowerWire® Pulse Radiofrequency Guidewire breaks or kinks.	Breaks and kinks in the PowerWire® Pulse Radiofrequency Guidewire are a potential cause of patient injury.	Discard immediately.

DISCLAIMER AND LIMITATION OF LIABILITY
THE LIMITED WARRANTY ABOVE IS THE SOLE WARRANTY PROVIDED BY SELLER. SELLER DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE OR PURPOSE.
THE REMEDY SET FORTH HEREIN SHALL BE THE EXCLUSIVE REMEDY FOR ANY WARRANTY CLAIM, AND ADDITIONAL DAMAGES, INCLUDING CONSEQUENTIAL DAMAGES OR DAMAGES FOR BUSINESS INTERRUPTION OR LOSS OF PROFIT, REVENUE, MATERIALS, ANTICIPATED SAVINGS, DATA, CONTRACT, GOODWILL OR THE LIKE (WHETHER DIRECT OR INDIRECT IN NATURE) OR FOR ANY OTHER FORM OF INCIDENTAL, OR INDIRECT DAMAGES OF ANY KIND, SHALL NOT BE AVAILABLE. SELLER'S MAXIMUM CUMULATIVE LIABILITY RELATIVE TO ALL OTHER CLAIMS AND LIABILITIES, INCLUDING OBLIGATIONS UNDER ANY INDEMNITY, WHETHER OR NOT INSURED, WILL NOT EXCEED THE COST OF THE PRODUCT(S) GIVING RISE TO THE CLAIM OR LIABILITY. SELLER DISCLAIMS ALL LIABILITY RELATIVE TO GRATUITOUS INFORMATION OR ASSISTANCE PROVIDED BY, BUT NOT REQUIRED OF SELLER HEREUNDER. ANY ACTION AGAINST SELLER MUST BE BROUGHT WITHIN EIGHTEEN (18) MONTHS AFTER THE CAUSE OF ACTION ACCRUES. THESE DISCLAIMERS AND LIMITATIONS OF LIABILITY WILL APPLY REGARDLESS OF ANY OTHER CONTRARY PROVISION HEREOF AND REGARDLESS OF THE FORM OF ACTION, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND FURTHER WILL EXTEND TO THE BENEFIT OF SELLER'S VENDORS, APPOINTED DISTRIBUTORS AND OTHER AUTHORIZED RESELLERS AS THIRD-PARTY BENEFICIARIES. EACH PROVISION HEREOF WHICH PROVIDES FOR A LIMITATION OF LIABILITY, DISCLAIMER OF WARRANTY OR CONDITION OR EXCLUSION OF DAMAGES IS SEVERABLE AND INDEPENDENT OF ANY OTHER PROVISION AND IS TO BE ENFORCED AS SUCH.
IN ANY CLAIM OR LAWSUIT FOR DAMAGES ARISING FROM ALLEGED BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, PRODUCT LIABILITY OR ANY OTHER LEGAL OR EQUITABLE THEORY, THE BUYER SPECIFICALLY AGREES THAT BMT SHALL NOT BE LIABLE FOR DAMAGES OR FOR LOSS OF PROFITS, WHETHER FROM BUYER OR BUYER'S CUSTOMERS. BMT'S LIABILITY SHALL BE LIMITED TO THE PURCHASE COST TO BUYER OF THE SPECIFIED GOODS SOLD BY BMT TO BUYER WHICH GIVE RISE TO THE CLAIM FOR LIABILITY.

XIV. LABELING AND SYMBOLS

	Manufacturer		Model number
	Use By		Do not reuse
	Caution		Lot Number
	Consult Instructions for Use		Non-Pyrogenic
	Sterilized using ethylene oxide		Keep Away from Sunlight
	Do Not Use if Packaging is Damaged		
	Follow Instructions for Use		
Rx ONLY	Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.		

XV. LIMITED WARRANTY – DISPOSABLES AND ACCESSORIES

Baylis Medical Technologies (BMT) warrants its Disposable and Accessory products against defects in materials and workmanship. BMT warrants that sterile products will remain sterile for a period of time as shown on the label as long as the original package remains intact. Under this Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMT will replace or repair, in its absolute and sole discretion, any such product, less any charges to BMT for transportation and labor costs incidental to inspection, removal or restocking of product. The length of the warranty is: (i) for the Disposable products, the shelf life of the product, and (ii) for the Accessory products, 90 days from shipment date.
This limited warranty applies only to new original factory delivered products that have been used for their normal and intended uses. BMT's Limited Warranty shall not apply to BMT products which have been re-sterilized, repaired, altered, or modified in any way and shall not apply to BMT products which have been improperly stored or improperly cleaned, installed, operated or maintained contrary to BMT's instructions.

No agent, employee or representative of Baylis Medical Technologies has the authority to bind the Company to any other warranty, affirmation or representation concerning the product.
This warranty is valid only to the original purchaser of Baylis Medical Technologies products directly from a Baylis Medical Technologies authorized agent. The original purchaser cannot transfer the warranty.
Use of any BMT product shall be deemed acceptance of the terms and conditions herein.

The warranty periods for Baylis Medical Technologies products are as follows:

Disposable Products	The shelf life of the product
Accessory Products	90 days from the shipment date

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