





Adelante[®] Breezeway[®] T

GUIDING SHEATH WITH HEMOSTATIC VALVE


The **Adelante® Breezeway®** guiding sheath is intended for the percutaneous introduction of various types of cardiovascular catheters to all chambers of the heart, including transseptal puncture. The Adelante® Breezeway® features a proprietary braided Flexsteer™ hydrophobic coated sheath and the latest SureSeal™ technology.




FEATURES

- Reliable Flexsteer™ reinforced braid technology provides exceptional control, kink resistance, while maintaining optimal flexibility
- Advanced SureSeal™ hemostatic valve technology to minimize air embolism and blood loss even under extreme conditions
- Proprietary hydrophobic coated sheath technology, enabling smooth access
- Tapered dilator to improve access
- Superior curve retention for reliable catheter placement
- Soft atraumatic distal tip with radiopaque marker band enabling accurate positioning under fluoroscopy
- Sideport holes on distal end for easy catheter aspiration and flushing
- Accepts maximum guidewire up to 0.035"/0.89 mm


SHAFT MODELS



SHAFT DIAMETERS



ORDERING INFORMATION



Set Includes:

- (1) Adelante Breezeway Sheath
- (1) Adelante Breezeway Dilator

Model # Set (One per box)	A1 Sheath Inner Diameter	A2 Usable Sheath Length	B1 Usable Dilator Length	C Compatible Transseptal Needle	Sheath Tip Curve
Straight Curve					
AB081210	● 8 F	61 cm	66 cm	71 cm	Straight
AB081211	● 8 F	79 cm	84 cm	89 cm	Straight
AB101218	● 10 F	61 cm	66 cm	71 cm	Straight
AB101219	● 10 F	79 cm	84 cm	89 cm	Straight
55° Curve					
AB081040	● 8 F	61 cm	66 cm	71 cm	55°
AB081044	● 8 F	79 cm	84 cm	89 cm	55°
AB101072	● 10 F	61 cm	66 cm	71 cm	55°
AB101076	● 10 F	79 cm	84 cm	89 cm	55°
70° Curve					
AB081041	● 8 F	61 cm	66 cm	71 cm	70°
AB081045	● 8 F	79 cm	84 cm	89 cm	70°
AB101073	● 10 F	61 cm	66 cm	71 cm	70°
AB101077	● 10 F	79 cm	84 cm	89 cm	70°
90° Curve					
AB081042	● 8 F	61 cm	66 cm	71 cm	90°
AB081046	● 8 F	79 cm	84 cm	89 cm	90°
AB101074	● 10 F	61 cm	66 cm	71 cm	90°
AB101078	● 10 F	79 cm	84 cm	89 cm	90°
120° Curve					
AB081043	● 8 F	61 cm	66 cm	71 cm	120°
AB081047	● 8 F	79 cm	84 cm	89 cm	120°
AB101075	● 10 F	61 cm	66 cm	71 cm	120°
AB101079	● 10 F	79 cm	84 cm	89 cm	120°

Guiding sheath supplied sterile. For single use only. Available bulk upon request. Standard lengths are noted above. Other lengths and curves available upon request.

Products not available in all regions.



P: 727.937.2511
F: 727. 934.9835
E: sales@oscor.com
www.oscor.com

CAUTION: Federal Law restricts this device to sale by or on the order of a physician (U.S.A.). Please refer to the Instructions for Use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.

