


Nykanen

RF Wire



A radiofrequency puncture wire used to cross delicate anatomy.

MAKE THE IMPACT OF A LIFETIME

Nykanen RF Wire




The **Nykanen RF Wire** is used to create controlled perforations in soft and delicate tissue.

CREATE CONTROLLED PUNCTURES WITH:



RF Energy

Deliver radiofrequency energy with the 0.016" active tip.




Flexible Wire

Flexible wire to pass through delicate anatomy.




Familiar Handling

Small 0.024" diameter, that handles like a guidewire.




Swift Exchange

Exchange quickly, with a free proximal end.



Atraumatic Tip

Atraumatic tip designed to avoid mechanical perforation.



Exchange Length

Ample wire length for exchange of catheters or tertiary devices.

SPECIFICATIONS

Feature	Specifications
Body Outer Diameter	0.024" (0.65mm)
Overall Length	265 cm
Active tip diameter	0.016"
Active tip length	1.5 mm

Recommended/compatible microcatheter - CIC38-145 (ProTrack™ Microcatheter)

Compatible Generator	Product Code
RFP-100A Generator*	RFK-265
RFP-100 Generator**	RFK-265

* Baylis Medical Company Radiofrequency Puncture Generator RFP-100A

** BMC Radiofrequency Puncture Generator RFP-100

U.S.A.: The Nykanen Radiofrequency Wire is intended to cut and/or coagulate soft tissues.
Canada: The Nykanen RF Wire is intended to create an atrial septal defect in the heart.

PRM-00304 EN J-1,2, V-5 © Copyright Baylis Medical Technologies Inc., 2007-2023. Baylis Medical Technologies Inc. reserves the right to change specifications or to incorporate design changes without notice and without incurring any obligation relating to equipment previously manufactured or delivered. The Baylis Medical Technologies logo is a trademark and/or registered trademark of Baylis Medical Technologies Inc. in the USA and/or other countries. Other trademarks are property of their respective owners. Patents Pending and/or issued. CAUTION: Federal Law (USA) restricts the use of these devices to or by the order of a physician. Before use, consult product labels and Instructions for Use for Indications for Use, Contraindications, Warnings, Precautions, Adverse Events and Directions for Use.

Products shown may not be approved in all jurisdictions.

Baylis

MED TECH

Baylis Medical Technologies Inc.
2645 Matheson Blvd. East
Mississauga ON Canada L4W 5S4

Tel.: 1 (888) 505-4885
www.baylismedtech.com / info@baylismedtech.com