

ORDERING GUIDE – Interventional Models

Navigate the vascular system with control and precision. The **Passport™** Steerable Sheath is the ultimate tool in gaining access to hard-to-reach target vessels.



Optimize your Workflow

Streamline Access

The deflectable tip facilitates access to challenging anatomies and minimizes the need for multiple fixed curve catheter exchanges.



Enhance Control

180° tip deflection enables smooth steering and navigation through the vascular system to target vessels.

Maximize Support

The reinforced braid design provides a stable yet flexible platform for delivery of diagnostic or therapeutic devices.


Multiple Diameter Options for Diverse Clinical Needs



ORDERING GUIDE – Interventional Models

Model # Set (One per box)	Sheath Inner Diameter	Usable Sheath Length	Usable Dilator Length	Curve Bend
STANDARD				
PSS0654509	●	6.5 F 45 cm	63 cm	9 mm
PSS0655507D2*	●	6.5 F 55cm	73 cm	7 mm
PSS0655509	●	6.5 F 55 cm	73 cm	9 mm
PSS0655517	●	6.5 F 55 cm	73 cm	17 mm
PSS0656509	●	6.5 F 65 cm	83 cm	9 mm
PSS0659009	●	6.5 F 90 cm	108 cm	9 mm
PSS0704509	●	7 F 45 cm	63 cm	9 mm
PSS0705509	●	7 F 55 cm	73 cm	9 mm
PSS0705517	●	7 F 55 cm	73 cm	17 mm
PSS0706509	●	7 F 65 cm	83 cm	9 mm
PSS0709009	●	7 F 90 cm	108 cm	9 mm
PSS0855517	●	8.5 F 55 cm	73 cm	17 mm
PSS1209030	●	12 F 90 cm	108 cm	30 mm
PSS1387139	●	13.8 F 71 cm	89 cm	39 mm
TD				
PSS0705509TD	●	7 F 55 cm	73 cm	9 mm

Shaft Configurations



Set Contains

- (1) Passport™ steerable guiding sheath
- (1) Passport™ straight curved dilator

* Contains second soft dilator with 80 cm usable length

