




ProTrack™
Microcatheter

FOR EXCELLENT TRACKABILITY.

ProTrack™ Microcatheter


Use to exchange guidewires, infuse diagnostic or therapeutic agents, and facilitate access to distal tortuous vasculature over a guidewire.

DESIGNED FOR VERSATILITY




Excellent Trackability

Facilitate access to distal, tortuous vasculature over a guidewire.



Easy Exchange

Exchange catheters effortlessly with a removeable hub.



Radiopaque Tip

Flexible 10 cm distal end with a radiopaque tip marker visible under fluoroscopy.

SPECIFICATIONS

Model	Outer Diameter Inches / French	Inner Diameter Inches	Length cm
CIC 35-145	0.035/2.7	0.022	145
CIC 38-145*	0.038/2.9	0.025	145

Maximum recommended pressure 250 psi

*Compatible with Nykanen RF Wire

PRM-00092 EN J-1,2 V-5 © Copyright Baylis Medical Technologies Inc., 2013-2023. Baylis Medical Technologies Inc. reserves the right to change specifications or to incorporate design changes without notice and without incurring any obligation relating to equipment previously manufactured or delivered. The Baylis Medical Technologies logo is a trademark and/or registered trademarks of Baylis Medical Technologies Inc. in the United States of America and/or other countries. Other trademarks are property of their respective owners. Before use, consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use. CAUTION: Federal Law (USA) restricts the use of this device to or by the order of a physician.

Products shown may not be approved in all jurisdictions.

Baylis MED TECH

Baylis Medical Technologies Inc.
2645 Matheson Blvd. East
Mississauga ON Canada L4W 5S4

Tel.: 1 (888) 505-4885
www.baylismedtech.com / info@baylismedtech.com