

ORDERING GUIDE – Interventional Models

Navigate the vascular system with control and precision. The **Passport®** Steerable Sheath is the ultimate tool in gaining access to hard-to-reach target vessels.



Optimize your Workflow

Streamline Access

The deflectable tip facilitates access to challenging anatomies and minimizes the need for multiple fixed curve catheter exchanges.


Enhance Control

180° tip deflection enables smooth steering and navigation through the vascular system to target vessels.

Maximize Support

The reinforced braid design provides a stable yet flexible platform for delivery of diagnostic or therapeutic devices.

Multiple Diameter Options for Diverse Clinical Needs




ORDERING GUIDE – Interventional Models

Model # Set (One per box)	Sheath Inner Diameter	Usable Sheath Length	Usable Dilator Length	Curve Bend
STANDARD				
PSS0654509	●	6.5 F 45 cm	63 cm	9 mm
PSS0655507D2*	●	6.5 F 55cm	73 cm	7 mm
PSS0655509	●	6.5 F 55 cm	73 cm	9 mm
PSS0655517	●	6.5 F 55 cm	73 cm	17 mm
PSS0656509	●	6.5 F 65 cm	83 cm	9 mm
PSS0659009	●	6.5 F 90 cm	108 cm	9 mm
PSS0704509	●	7 F 45 cm	63 cm	9 mm
PSS0705509	●	7 F 55 cm	73 cm	9 mm
PSS0705517	●	7 F 55 cm	73 cm	17 mm
PSS0706509	●	7 F 65 cm	83 cm	9 mm
PSS0709009	●	7 F 90 cm	108 cm	9 mm
PSS0855517	●	8.5 F 55 cm	73 cm	17 mm
PSS1209030	●	12 F 90 cm	108 cm	30 mm
PSS1387139	●	13.8 F 71 cm	89 cm	39 mm
TD				
PSS0705509TD	●	7 F 55 cm	73 cm	9 mm

Shaft Configurations

Pre-formed Shaft (TD)



Straight Shaft (Standard)



Set Contains

- (1) Passport® steerable guiding sheath
- (1) Passport® straight/curved dilator

* Contains second soft dilator with 80 cm usable length



DISTRIBUTED BY –
BAYLIS MEDICAL TECHNOLOGIES

2645 Matheson Blvd. East
Mississauga, ON, CA
L4W 5S4

Work: 905-948-5800 x334
Toll Free: 1-888-505-4885

www.baylismedtech.com



MANUFACTURED BY –
OSCOR INC.

3816 DeSoto Blvd.
Palm Harbor, FL 34683 U.S.A

www.oscor.com



PRM-00882 EN J-1, 2, V-2 © 2024 Baylis Medical Technologies Inc. Passport and the Baylis Medical Technologies Logo are trademarks and/or registered trademarks of Baylis Medical Technologies Inc. **CAUTION:** Federal Law (USA) restricts the use of the devices to or by the order of a physician. Before use, consult product labels and Instructions for Use for Indications for Use, Contraindications, Warnings, Precautions, Adverse Events and Directions for Use.

Products shown may not be approved or available in all jurisdictions.