

Passport[®]

Steerable Sheath

Large Bore

Large bore access with
precise steering




For product
information and
clinical references,
visit our website.

Baylis
MEDTECH

Advancing Technology.
Transforming Care.™

Large Bore Sheaths

Large bore sheaths provide increased inner diameter to accommodate larger devices or multiple tools with the enhanced control and stability of steerable technology.



Product Code	Sheath Inner Diameter	Usable Sheath Length	Usable Dilator Length	Curve Bend
--------------	-----------------------	----------------------	-----------------------	------------

PASSPORT® STEERABLE SHEATHS

PSS1209030	● 12 F	90 cm	108 cm	30 mm
PSS1387139	● 13.8 F	71 cm	89 cm	39 mm

DESTINO™ REACH STEERABLE SHEATHS

DCR1207522	● 12 F	75 cm	93 cm	22 mm
DCR1207150	● 12 F	71 cm	89 cm	50 mm
DCR1207122	● 12 F	71 cm	89 cm	22 mm

DESTINO™ STEERABLE SHEATH SETS

DS126722	● 12 F	67cm	90 cm	22 mm
----------	--------	------	-------	-------

DESTINO™ STEERABLE SHEATH KITS

DS126750GW	● 12 F	67 cm	90 cm	50 mm
DS127122GW	● 12 F	71 cm	94 cm	22 mm

Passport® steerable sheaths include: (1) Passport® steerable guiding sheath, (1) Passport® straight dilator

Destino™ Reach steerable sheaths include: (1) Destino™ Reach steerable guiding sheath, (1) 50° curved dilator

Destino™ steerable sheaths sets include: (1) Destino™ steerable guiding sheath, (1) Straight dilator

Destino™ steerable sheaths kits include: (1) Destino™ steerable guiding sheath, (1) Straight dilator, (1) Guidewire outer diameter: 0.035", L: 180 cm